Instructions for Use EN

■ PRODUCT DESCRIPTION
Au-based dental casting alloy, Typ 4

■ INDICATIONS *
Inlays, Onlays, Crowns, Telescopic Crowns, Conus Crowns, Bridges, Wide Bridges, Cast Posts / Cores, Bars, Attachments, Inplant Retained Superstructures, Partial Dentures

■ WAXING / MODELLATION
Design the framework in full contour or in a reduced anatomic shape taking the planned veneer into consideration. For Lab Composite veneering material use mechanical retentions. Single crowns require a minimum thickness of 0.3 mm. Abutment crowns require a minimum thickness of 0.5 mm. Ensure the specific technique used for the veneering material. Avoid sharp angles. Connectors must have the required dimensions to provide resistance to deformation. Create large surface areas for planned soldering, with a gap of 0.05 – 0.2 mm.

■ SPRUING
Provide the modeled single-tooth restoration or bridge framework with sprues of a suitable size. In general the reservoir, sprue leads, and connector sprues, whether pear shaped or traditional, must be sized according to the specific technique used. When using the direct or indirect technique be sure that the reservoir is positioned in the heat center. The connector sprues between the reservoir and the casting should be a maximum of 2.5–3.0 mm in length and width. The wax pattern including the sprues must be weighed in grams in order to determine the needed amount of alloy. Wax conversion formula: wax weight (gram) x alloy density = grams of alloy required.

■ INVESTING
Use a phosphate- / gypsum-bonded investment material. Follow the manufacturer’s instructions.

■ PREHEATING / BURN-OUT
Recommended burn-out temperature: 700 °C

■ MELTING AND CASTING
Torch: Propane 0.15 bar; Oxygen 0.35 bar

Other specifics may be required by the type of casting machine. It is recommended to use a separate and clean carbon /ceramic crucible for each alloy. The recommended ratio of used material to new material is 1:1. Use casting flux if needed.

■ Casting Temperature: 1060 °C

■ FRAMEWORK FINISHING
After benching, carefully divest and clean the casting with aluminum oxide (AL₂O₃). Do not use a hammer for divesting. Finish the casting with carbide burs and/or with ceramic-bonded grinding instruments. For veneering with Lab Composite material, the framework must be finished and polished. Subsequently, steam clean or ultrasonic clean with distilled water or ethanol and dry the framework. Use the appropriate Lab Composite material following the manufacturer’s instructions.

■ HEAT TREATMENT
Softening: 600 °C for 15 min; bench cool
Hardening: 400 °C for 15 min; bench cool

■ SOLDERING AND LASER WELDING

The soldering gap should not be wider than the thickness of the soldering material. Allow the soldered casting to cool slowly. Use flux sparingly.

Pre Solder:	Auropal Lot M-1	Flux:	Bondal Flux
Post Solder:	Auropal Lot W-2	Flux:	Bondal Flux
Laser Welding Wire:	Portadur P2		

■ POLISHING
Carefully remove any oxide and flux residue. Smooth the metal surfaces with rubber polishes. Polish to a high gloss finish using polishing paste. Subsequently, clean using ultrasonic cleaning equipment or careful steam cleaning.

ADDITIONAL SAFETY CONCERNS AND INSTRUCTIONS

■ CONTRAINDICATIONS
Patients with allergy/sensitivity to any major or minor elements of this alloy, consultation with a physician is recommended. Alloy is not to be used for any application not included within the indications.

■ SIDE EFFECTS

In individual cases, sensitivity or allergies to elements of this alloy may occur. Ivoclar Vivadent makes no claims regarding the MRI-compatibility of its dental alloys. It is recommended that the patient be made aware of the possibility for dental alloys to affect MRI results and to disclose the presence of dental alloys to the MRI technician prior to conducting a test.

■ INTERACTIONS

Galvanic effects may occur between different or dissimilar alloys in the same oral environment.

■ CAUTION
Metal vapors and metal dust are harmful if inhaled. Therefore, the use of extraction equipment and/or a suitable protective mask is advised!

■ STORAGE CONDITION
Store in a dry environment at room temperature.

■ DISCLAIMER

This material has been developed solely for use in dentistry. Processing should be carried out strictly according to the instructions. The user assumes full responsibility for any damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the products for their suitability and use for any purpose not explicitly stated in the Instructions. These regulations also apply if the materials are used in conjunction with products of other manufacturers.

PROCESSING DATA			
Investment Material:	phosphate-bonded / gypsum-bonded		
Preheating / Burn-out Temperature:	700 °C		
Crucible:	Carbon crucible / Ceramic crucible		
Casting Temperature:	1060 °C		
Recommended Lab-Composite:	SR Nexos®		
Softening:	600 °C for 15 min; bench cool		
Hardening:	400 °C for 15 min; bench cool		
Pre Solder / Flux:	Auropal Lot M-1	Bondal Flux	
Post Solder / Flux:	Auropal Lot W-2	Bondal Flux	
Laser Welding Wire:	Portadur P2		
TECHNICAL DATA (ISO 22674:2016)			
Type / Color:	4	Rich Yellow	
Density (g/cm ³):	14,3		
Melting Range (Solidus/Liquidus):	870–910 °C		
Elastic Modulus (GPa):	101		
	As Cast	Softened	Hardened
Vickers-Hardness:	260	160	260
Tensile Strength (MPa):	740	480	740
0,2% Proof Stress (MPa):	690	350	680
Elongation (%):	9	34	10
* Siehe TYP-KLASSIFIZIERUNG DUE TO PHYSICAL PROPERTIES			

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