## Instructions for Use

**■ PRODUCT DESCRIPTION**  
Ag-based dental casting /metal-ceramic alloy, Type 4

**■ INDICATIONS\***  
Inlay, Onlays, 3/4 Crowns, Crowns, Telescope Crowns, Conus Crowns, Bridges, Cast Posts/Cores, Bars, Attachments, Implant-Suprastructures, Partial Dentures

**■ WAXING/MODELATION**

Design the framework in full contour or in a reduced anatomy shape taking the planned veneer into consideration. Use the compatible waxing material as recommended. Single crowns require a minimum thickness of 0.3 mm. Abutment crowns require a minimum thickness of 0.5 mm. Ensure the correct design provides adequate support for the veneering material. Avoid sharp angles. Connectors must have the required dimensions to provide resistance to deformation. Create large surface areas for planned soldering, with a gap of 0.05–0.2 mm.

**■ SPRUING**  
Provide the modelled single-tooth restoration or bridge framework with sprues of a suitable size. In general the reservoir, sprue leads, and connector sprues, whether pear shaped or traditional, must be sized according to the spruing technique used. When using the direct or indirect technique be sure that the reservoir is positioned so that the thickness of the wax between the reservoir and the casting should be a maximum of 2.5–3.0 mm in length and width. The wax pattern including the sprues must be weighed in grams in order to determine the needed amount of alloy. Wax conversion formulae: wax weight (grams) / alloy density = grams of alloy required

**■ INVESTING**  
Use a phosphate-bonded investment material. Follow the manufacturer's instructions.

**■ PREHEATING / BURN-OUT**  
Recommended burn-out temperature: 800 °C

**■ MELTING AND CASTING**  
Torch: Propano 0.15 bar, Oxygen 0.35 bar  
Note: The correct wax may be required by the type of casting machine. It is recommended to use a separate and clean carbon/ceramic crucible for each alloy. Preheat the ceramic crucible in the burnout furnace. Do not use used material to new material 1:1. Do not use flux.

**■ HEAT TREATMENT**  
Casting Temperature: 1320 °C  
**■ FRAMEWORK FINISHING**  
After being cooled, carefully invest and cast the casting with aluminum oxide (Al<sub>2</sub>O<sub>3</sub>). Do not use a hammer or other specific tool which may be required by the type of casting machine. It is recommended to use a separate and clean carbon/ceramic crucible for each alloy. Preheat the ceramic crucible in the burnout furnace. Do not use used material to new material 1:1. Do not use flux.

**■ POLISHING**  
Remove any oxide and flux residue. Smooth the metal surfaces with rubber polishers. Polish to a high gloss finish using polishing paste. Subsequently, clean using ultrasonic cleaning equipment or careful steam cleaning.

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### ADDITIONAL SAFETY CONCERNS AND INSTRUCTIONS

**■ CONTRAINDICATIONS**  
For patients with known allergy/sensitivity to any major or minor elements of this alloy, consultation with a physician is recommended. Alloy is not to be used for any application not included within the indications.

**■ SIDE EFFECTS**  
In some individuals, sensitivity or allergies to elements of this alloy may occur. Ivoclar Vivadent makes no claims regarding the MRI-compatibility of its dental alloys. It is recommended that the patient be made aware of the possibility of facial distortion as a result of MRI results and to disclose the presence of dental alloys to the MRI technologist prior to conducting a test.

**■ INTERACTIONS**  
Galvanic effects may occur between different or dissimilar alloys in the same oral environment.

**■ CAUTION**  
Metal vapors and metal dust are harmful if inhaled. Therefore, the use of extraction equipment and /or suitable protective masks is advised

**■ STORAGE CONDITION**  
Store in a dry environment at room temperature.

**■ DISCLAIMER**  
This material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the products for their suitability and for any purpose not explicitly stated in the Instructions. These regulations also apply if the materials are used in conjunction with products of other manufacturers.

PROCESSING DATA			
Investment material:	phosphate-bonded		
Preheating / Burn-out Temperature:	800 °C		
Crucible:	Carbon crucible / Ceramic crucible		
Casting Temperature:	1320 °C		
Oxidation:	Temperature: 800 °C; Holding time: 5 min; Vacuum: No		
CTE:	(25–500 °C): 16,7 x 10 <sup>-6</sup> /K		
Recommended Veneering Materials:	low-fusing, high-expanding, leucite-containing metal ceramics		
Hardening:	600 °C for 15 min; bench cool		
Pre Solder / Flux:	Porta IP Lot V-1 High Fusing Bondal Flux		
Post Solder / Flux:	Porta OP Lot W-2 Bondal Flux		
Laser Welding Wire:	Porta P6		

TECHNISCHE DATEN (ISO 22674:2016 & ISO 9693-1:2012)			
Typ / Color:	4 White		
Density (g/cm <sup>3</sup> ):	10,7		
Schmelzbereich (Solidus/Liquidus):	1100–1170 °C		
Elastic Modulus (GPa):	104		
As Cast	Porcelain Fired	Hardened	
Viickers-Härde:	190	210	210
Tickers-Härde (MPa):	190	210	210
Resistenz alla trazione (MPa):	380	380	380
0,2% Dehnungsreaktion (MPa):	380	380	380
Bruchdehnung (%):	25	23	22

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